



## **MULTI-FREQUENCY VIBRATOR**

VESTIBULAR VIBRATOR

### **TECHNICAL SPECIFICATIONS**

- **Vibration frequency:** 30, 60 or 100 Hz  $\pm$  5%
- **Vibration amplitude:** 1 mm
- **Power supply:** Control box 24 V powered
- **Maximum power:** 20 VA
- **External power supply:** 24 V – 1,25 mA
- **Medical device (Directive 93/42/EEC):** Class IIa
- **Standard:** EN 60601-1 / EN 60601-1-2

### **INCLUDED PARTS**

- Control box x 1
- External power supply x 1
- Biocompatible foam cover x 10

### **SPARE PARTS**

- Biocompatible foam cover x 30

### **PRODUCT CODES**

**10741:** Multi-frequency vibrator (30, 60, 100 Hz) - SY.VIBREUR\_II

### **MANUFACTURER**

**Synapsys S. A.**

2, rue Marc Donadille  
13013 Marseilles (France)

### **INVENTIS S.R.L.**

CORSO STATI UNITI, 1/3 | TEL: +39 049.8962 844  
35127 PADOVA - ITALIA | FAX: +39 049.8966 343

C.F./P.IVA:IT03957810280 |  
R.E.A. PD 350789

INFO@INVENTIS.IT  
WWW.INVENTIS.IT