



DIAGNOSIS

MED4 chair

Electronic rotary chair

Synapsys has designed the MED4 chair for an ENT practice as to complete the VNG tests with the kinetic rotational test.

KEY BENEFITS

It enables to calculate:

- the vestibular gain
- the phase
- the symmetry
- the time constants

It offers a large range of kinetic tests:

- Various sinusoidal tests (SHA profiles, decreasing sinus)
- Multi-frequencies (sweep)
- Multi-turn (step velocity)
- Sudden stop stimulations (square back and forth)

Its ergonomics ensures the **patient's comfort and security** when carrying out middle frequencies tests.

The MED4 enables the patient to be in different positions and thus can be used for several vestibular tests complementary to the kinetic test.



OPTIONS:

Safety harness

RELATED PRODUCT:

VNG Ulmer

CONNECTIVITY:

The MED4 can be controlled by VNG software or used manually.

