



# Maestro

SOFTWARE

## DESCRIPTION

Maestro is a software for managing patient data that can interface with all the Inventis device, controlling or getting exams results from them.

## OPERATING MODES

- Noah module
- Stand-Alone with proprietary database
- Stand-Alone without database functions (file based)
- Third-party Applications (exchange data with external applications): MII proprietary protocol and GDT standard

## AVAILABLE TESTS

- Audiology exams
- REM
- HIT
- Video-Otoscopy
- Video and Images acquisition
- Questionnaire
- Remote audiometry and video-otoscopy

## COMPATIBLE DEVICES

- Triangle (Audiometer)
- Bell (Audiometer)
- Harp (Audiometer)
- Piano (Audiometer)
- Piccolo (Audiometer)
- Cello (Audiometer)
- Trumpet (REM, Audiometer)
- Drum (HIT)
- Flute (Middle-Ear Analyzer)
- Clarinet (Middle-Ear Analyzer)
- Viola (Middle-Ear Analyzer, Audiometer)
- Timpani (Handheld Tympanometer, Audiometer)
- Delfino (Otoscope)
- MediCam (Video-endoscopy camera)
- Analog-to-USB video converter

## INTERFACE WITH EQUIPMENT

USB connection

## PC MINIMUM REQUIREMENTS

CPU: Intel® i5, 6<sup>th</sup> generation or above, 8 GB RAM

Display:

- Resolution: 1366x768 (16:9) or higher
- Scale and Layout: 100% dpi (recommended)\*

Operating System: Windows 10 and Windows 11, 32 or 64 bits

\*: higher scale settings can be set, but not with the lowest display resolutions

## NOAH MINIMUM REQUIREMENTS

NOAH 4.4 or greater (or NOAH Engines supporting NOAH 4.4 or greater)

## SOFTWARE SUPPORT

Remote assistance tool: included

## SOFTWARE UPDATES

Available for free

## MII INTERFACE

Library for interfacing Maestro with 3<sup>rd</sup> party software

Available upon request (IT technician intervention requested)

Maestro is developed by:

### INVENTIS s.r.l.

CORSO STATI UNITI, 1/3

35127 PADOVA – ITALY

PHONE: +39.049.8962 844

FAX: +39.049.8966 343

info@inventis.it

www.inventis.it

Inventis® is a registered trademark of Inventis s.r.l.