

FITTING

# Video Otoscope H

Video Otoscope



  
**Orchestral perfection**  
A full line of Audiology Equipment

- High Quality Images (Full-HD)
- Built-in display
- Live view mode with Maestro
- Suitable for remote video-otoscopy

# Video Otoscope H

Video Otoscope



Video Otoscope H is a digital handheld video otoscope, ergonomic and compact and extremely easy to carry. The hand piece integrates the camera, the LED illumination, the display, and the rechargeable battery. It is fully integrated with the software Maestro: it can be used in connection with the PC to store the images directly in the database, or as a standalone device and download the acquired images later. The live image can be seen both in Maestro and in the integrated display. It can also be used with the software Satellite for remote video otoscopy.

## TECHNICAL DATA

### VIDEO OUTPUTS CHARACTERISTICS

Digital: USB 2.0 (mini-USB)  
Composite AV

### USE WITH MAESTRO AND SATELLITE

- Live image and acquisition
- Data transfer
- Remote video otoscopy

### POWER SUPPLY

Output: 5V DC, 1.2A  
Input: 100-240 V AC, 50/60 Hz

### MECHANICS

Size (WxHxD): 7 x 18 x 7 cm / 2.8 x 7.1 x 2.8 in  
Weight: 200 g / 7.1 oz

### APPLICABLE STANDARDS

Electrical safety: EN 60601-1  
EMC: EN 60601-1-2

### MDR CLASSIFICATION

MDR (EU) 2017/745 classification: Class I (Rule 13)

### MINIMUM COMPUTER REQUIREMENTS

Refer to Maestro software technical specifications for indications about the minimum system requirements.

 Find our distributors at:  
[www.inventis.it/dealers](http://www.inventis.it/dealers)

Discover the complete orchestra at: [www.inventis.it](http://www.inventis.it)

