

All you need for the perfect fitting



A complete line of innovative and unique solutions



Reliability and ease of use



REM system, Video otoscopy, HIT box, counseling software

TRUMPET

REM SYSTEM-AUDIOMETER



Trumpet is an advanced and easy-to-use Real Ear Measurement system, capable of performing quick and accurate in-situ measures thanks to its small and lightweight probes and its built-in amplified speaker. Not only, Trumpet combines also a diagnostic audiometer performing pure-tone and speech audiometry.

- Standard meas. (REUG/R; REAG/R; REIG; REOG/R)
- Speech mapping
- PTA and Speech audiometry, both in AC, BC, and FF (AUD version)
- Wireless probes (opt.)

HARMONICA

VIDEO-OTOSCOPE



A handy digital video-otoscope with chip-on-tip high definition camera. Fully integrated with Maestro software and compatible with Noah, **Harmonica** video-otoscope offers the possibility to capture outstanding quality images of the eardrum, that can be shown in real-time to the patient as well.

- Chip-on-tip HD camera (resolution: 1280 x 720)
- Crisp and clear images
- Noah compatible



The integration between the diagnostic devices and the software is essential to optimize the daily practice while maintaining the concentration on the patient. With **Maestro**, a few simple clicks are enough to store all the instrumental exams in a single folder, also accessible remotely. Maestro is compatible with Noah 4.

Maestro



DRUM

HIT BOX

Drum is a hearing instrument test box (HIT) capable of an impressive ambient sounds attenuation, and featured with all the standard tests for evaluating performances of air conduction hearing aids.

- High attenuation of ambient sounds
- Standard tests according to IEC 60118-7 and ANSI S3.22
- Battery adapters
- Built-in telecoil

SYMPHONIA

VIRTUAL SOUND ENVIRONMENT CREATION



Symphonia allows to simulate a realistic sound environment and makes the patient feel the benefits provided by advanced hearing aid features, without having to leave the room. Symphonia's advanced algorithm reproduces background noises recorded with special multichannel microphones, to make the patient feel inside a realistic external environment.

- Realistic sound environment
- Directional sources
- Flexible loudspeakers configuration
- Advanced counseling tool

Fitting

Including video-otoscopes, a real-ear measurement system, a hearing instrument test box and a software to simulate real-life sound environments, this product range is all you need for a perfect fitting.

Fitting



Trumpet

REM System-Audiometer

Fitting



Harmonica

Video otoscope

Fitting



Drum

HIT Box

Fitting



Symphonia

Virtual sound environment creation

Find your distributor:

www.inventis.it



academia

Webinar. Knowledge.Inventis.

Follow us on social media to stay updated on all academia appointments.

ITALY

INVENTIS S.R.L. A SOCIO UNICO
CORSO STATI UNITI, 1/3
35127 PADOVA
TEL.: +39.049.8962.844

NORTH AMERICA

INVENTIS NORTH AMERICA INC.
2503 S WASHINGTON AVE #586
TITUSVILLE FL, 32780
TOLL FREE: +1.844.683.6847

FRANCE

SYNAPSYS SAS
2 RUE MARC DONADILLE
HÔTEL TECHNOPTIC
13013 MARSEILLE
TÉL.: +33.4.91.11.75.75